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
**SUPERSECT/ SUPERLOOP PK FRONT
LOADING ELECTRODES**

RESECTOSCOPE INSTRUMENTS

REF. 784515 SuperSect PK® FRONT LOADING ELECTRODE

REF. 784415 SuperLoop PK® FRONT LOADING ELECTRODE

IMPORTANT: These Instructions For Use should be retained for reference until the contents of the multi-pack product package are exhausted. They should be used in conjunction with the User Manual provided with the Gyrus ACMI PlasmaKinetic SuperPulse Generator.

Instructions for use

IMPORTANT: Refer to the relevant product Instructions For Use and the User Manual for full details of system preparation other than for information relating directly to the PlasmaKinetic SuperSect/Loop device.

Indications

The Gyrus ACMI PlasmaKinetic SuperSect/Loop devices are bipolar instruments intended for use in urological electrosurgical procedures involving the ablation or removal of soft tissue and where associated haemostasis is required. The specific urological indications are:

- Transurethral Prostatectomy (TURP) for benign prostatic hypertrophy,
- Resection of bladder tumours,
- Cystodiathermy.

Before Surgery

1. Introduce a new PlasmaKinetic SuperSect/Loop device and the instrument cable into the sterile field. Insert the T-shaped connector on the instrument cable into the slot on the base of the white instrument mounting block. Ensure that the connector is fully seated within the block.

Insert a telescope into the working element. Load the instrument from the front (Fig 1), inserting the proximal end into the working element past the locking mechanism. Attach the stabilizer clip to the underside of the telescope. Gently pull the stabilizer clip to verify that the electrode proximal end is securely fastened.

2. Continue to assemble the resectoscope according to the details provided in the User Manual.
3. Connect the Instrument to the generator by inserting the 3-pin generator connector into the PK/SP outlet on the generator console. Once the connection is made, the "CONNECT PK CABLE" message on the Generator display will change to the default settings appropriate to the instrument classification code.

CAUTION: Verify that the Gyrus ACMI PlasmaKinetic SuperPulse Generator display matches the PlasmaKinetic SuperSect/Loop device default settings on the relevant package labelling. If the Generator display does not match the default settings then please contact Gyrus ACMI or your local authorised service representative.

Visually inspect the outer surfaces of the resectoscope and device for unintended roughness, sharp edges or protrusions that may be the cause of a safety hazard.

CAUTION: Failure to follow the correct assembly sequence could result in an inadvertent application of excess RF energy.

During Surgery

Refer to the User Manual for instructions concerning surgical use of the electrode.

After Surgery

1. Detach the instrument from the generator by removing the 3-pin generator connector from the console.

CAUTION: If changing Electrodes during the procedure always use a fresh cable that is supplied with the Electrode; failure to use the fresh cable supplied with the Electrode can result in misidentification of the Electrode by the generator.

CAUTION: Always detach by grasping the connector and pulling outward and not the cable. Failure to do so may result in product damage.

2. Disassemble the sheaths from the working element. Remove the instrument from the working element by depressing the locking mechanism button on the white instrument block. With the button fully depressed, remove the instrument from the working element. Remove the T shaped connector from the instrument block on the working element.

CAUTION: Only remove the instrument when the button is fully depressed. Removal of the instrument without depressing the button fully can result in damage to the working element.

3. It is generally advisable to inspect electrosurgical instruments for integrity after use.
4. Dispose of the SINGLE USE instrument(s) and Instrument cable(s) according to your facility's policy on the destruction of surgical waste.

WARNING: Do not attempt to re-use PlasmaKinetic Instruments or Instrument cable. Heat or chemical sterilisation may render the instrument or Instrument cable mechanically or electrically unsafe.

GENERAL WARNINGS AND PRECAUTIONS

Contraindications for Use

The use of this device is contraindicated in patients with the following conditions:

- Carcinoma of the bladder or prostate without tissue diagnosis.
- Patients with urinary tract infection.
- Patients with incipient renal failure.

It is advisable to monitor the input and output volumes of the irrigation fluid in all patients but especially those with cardiovascular insufficiency or poor renal function.

Patients with Pacemakers

Use with caution in the presence of internal or external pacemakers. Interference from an electrosurgical current can cause a pacemaker to enter an asynchronous mode or can block the pacemaker effect entirely. For further information, consult the pacemaker manufacturer or hospital Cardiology Department.

IMPORTANT: Sterility is guaranteed unless the package is opened or damaged. Do not resterilize. If product is opened or damaged it should be discarded and not used.

Part No. 184007-EE

- 1 -

CAUTION: Use only Gyrus ACMI approved accessories supplied for use with this product. Product damage or accessory failure may otherwise result during use.

WARNING: The use of non Gyrus ACMI Electrodes with the Gyrus ACMI instrument cable is not permitted, use of non supported Electrodes may result in user or patient harm.

WARNING: Use extreme caution when using electrosurgery in close proximity to or in direct contact with any metal objects. The working channel and operating sheaths of most rigid endoscopes are metal. Do not activate the instrument while any portion of the instrument tip is within the working channel or in contact with another metal object. Localised heating of the instrument and the adjacent metal object or working channel may result in damage to the contacting endoscope, and/or instrument tip.

WARNING: When not in use, place active instrument in a clean, dry, non-conductive, and highly visible area not in contact with the patient. Inadvertent activation while in contact with the patient may result in burns.

WARNING: Should the device become deformed or damaged in anyway, it should be replaced immediately as continued use may result in unpredictable operation or device failure.

WARNING: Fire Hazard: Do not place active accessories near or in contact with flammable materials (such as gauze or surgical drapes). Electrosurgical accessories which are activated or hot from use can cause a fire.

WARNING: The PlasmaKinetic Instrument must be used with saline as an irrigation fluid. Saline is critical to the proper functioning of the instrument. The use of other irrigation fluids will result in nonperformance of the instrument and may result in injury to the patient.

WARNING: Risk of injury to patient: Thermal cell damage can occur when using preheated saline. Always make sure that the temperature of the saline does not rise above body temperature (37 °C/99 °F).

WARNING: Risk of injury to patient: Localized excessive heating of the saline can be caused by HF current. Always make sure to use a sufficient flow rate (minimum 1 litre/5 min).


CAUTION: Proper use of the system relies on tissue contact with the active tip of the instrument only. Do not bury the tip in tissue beyond the junction of the active instrument and the ceramic insulator as product damage may result during use.

This product is for single-use only.

It has not been designed to be re-used or re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion and dulled edges, ceramic and plastic deformation or splitting which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.

Key to symbols used on package labelling

	- Sterilized by radiation		- Storage Temperature Range (°C)
	- Lot Number		- Single use only
	- Keep away from Direct Sunlight		- Use by date shown
	- Keep Dry		
	- Fragile, Handle with Care		- When a product is manufactured for Gyrus ACMI
	- Represents the quantity of saleable units inside the package.		
	- Represents the quantity of devices inside the package.		
	# - Numeral corresponds to the number of devices inside the		

Follow the correct sequence of assembly and disassembly of the cable and electrode to avoid damage

Manufactured in the UK

Manufactured For:
Gyrus ACMI Inc.
136 Turnpike Road
Southborough, MA 01772-2104, USA


Customer Service: 1 800 852 9361
Technical Service: 1 800 621 3739

EC REP

Gyrus Medical Limited,
Fortran Road, St Mellons
Cardiff CF3 0LT, UK.

PlasmaKinetic™, PK® and SuperPulse® are trademark(s) of Gyrus Medical, Ltd. and/or its affiliated entities, in the U.S., and/or other countries.

Rx only - CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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